

Technical Data Sheet

Product Name: TGF b 2 Recombinant Protein

Catalogue Number: BPBOI052

Volume: 10µg-10mg

Description:

Source : Nicotiana benthamiana. TGFB2 Human Recombinant produced in plants is a homodimeric polypeptide chain containing 2 x 118 amino acids and having a total molecular mass of 27.08kDa. The TGFB2 is fused to 6xHis Tag at Nterminus and purified by proprietary chromatographic techniques. TGFB2 is a 27.08 kDa protein having two identical 118 amino acid peptide chains linked by a single disulfide bond. TGFB2 is part of a family of five related cytokines that have an extensive variation of normal and neoplastic cells, indicating the importance of these homo-dimer proteins as multifunctional regulators of cellular activity. The three mammalian isoforms of TGF-? (TGFb1, TGFb2 and TGFb3) signal through the same receptor and stimulate similar biological responses. They are involved in physiological processes as embryogenesis, tissue remodelling and wound healing.

Product Info:

Amount : 5 µg

Purification : Greater than 97.0% as determined by SDS-PAGE.

Content : Lyophilized from a concentrated (1mg/ml) solution containing 50mM Tris-HCl pH-7.4.

Storage condition : Lyophilized TGFB2 although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution TGFB2 Human should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

Amino Acid: HHHHHHALDAAYCFRNVQDNCCLRPLYIDFKRDL
GWKWIHEPKGYNANFCAGACPYLWSSDTQHSLVLSLYN
TINPEASASPCCVSQDLEPLTI
LYYIGKTPKIEQLSNMIVKSKCS.

Application Note

It is recommended to reconstitute the lyophilized TGFB2 in sterile 18M-cm H₂O not less than 1Åµg/40Åµl, which can then be further diluted to other aqueous solutions. The biological activity of TGFB2 is measured in culture by its ability to inhibit the mink lung epithelial (Mv1Lu) cells proliferation. ED₅₀ < 40ng/ml, corresponding to a specific activity of 25,000 units/mg.



An ISO 9001:2015 & GMP Certified Company

Technical Assistance:

At BioPioneer, we are committed to providing top-notch technical support and ensuring its accessibility. For any technical assistance, please email us at info@biopioneer.in.