

Technical Data Sheet

Product Name: PLGF 1 Recombinant Protein

Catalogue Number: BPBOI098

Volume: 10µg-10mg

Description:

Source : Insect Cells. Placenta Growth Factor-1 Human Recombinant produced in insect cells is a homodimer, glycosylated polypeptide chain containing 2 x 131 amino acids and having a total molecular mass of approximately 34 kDa. The PLGF-1 is purified by proprietary chromatographic techniques. PLGF is a growth factor active in angiogenesis, and endothelial cell growth, stimulating their proliferation and migration. It binds to receptor vegfr-1/flt1.

Product Info:

Amount : 10 µg

Purification : Greater than 95.0% as determined by SDS-PAGE.

Content : Lyophilized from a concentrated (1mg/ml) solution containing 50mM acetic acid.

Storage condition : Lyophilized Placenta Growth Factor 1 although stable at room temperature for 3 weeks, should be stored desiccated below – 18°C. Upon reconstitution PLGF1 should be stored at 4°C between 2-7 days and for future use below -18°C. Please prevent freeze-thaw cycles.

Amino Acid : The sequence of the first five N-terminal amino acids was determined and was found to be Met- Leu-Pro-Ala-Val.

Application Note

It is recommended to reconstitute the lyophilized Placenta Growth Factor 1 in sterile 0.1M acetic acid not less than 100Å, Åµg/ml, which can then be further diluted to other aqueous solutions. The ED50 as determined by its ability to bind to immobilized rh-sFlt-1 in a functional ELISA. PLGF-1 human Recombinant can bind to immobilized rh-sFlt-1 (100ng/well) with a linear range at 0.5-10ng/ml, corresponding to a Specific Activity of 1x10⁵-2x10⁶units/mg.

Technical Assistance:

At BioPioneer, we are committed to providing top-notch technical support and ensuring its accessibility. For any technical assistance, please email us at info@biopioneer.in.

